PLAIN / CHROMIC GUT

Instructions for Use

PLAIN AND CHROMIC GUT (ABSORBABLE SURGICAL SUTURES, USP / EP)

DESCRIPTION

Plain and Chromic Gut is an absorbable sterile surgical suture composed of purified connective tissue (mostly collagen) derived from either the serosal layer of beef (bovine) or the submucosal fibrous layer of sheep (ovine) intestines. The sutures are available plain or chromic and packaged in a 89% isopropanol, 10% water, and 1% triethanolamine tubing fluid. The Chromic Gut has been tanned by treatment with chromium salts in trivalent form and with oxidised pyrogallol (oxidised pyrogallic acid) to prolong its resistance to absorption.

Plain and Chromic Gut meets all requirements established by the United States Pharmacopeia (USP) and European Pharmacopeia (E.P) for Absorbable Surgical Sutures.

INDICATIONS

Plain and Chromic Gut is indicated for use in general soft tissue approximation and/or ligation, including use of ophthalmic procedures, but not for use in cardiovascular and neurological procedures.

ACTIONS

When Plain and Chromic surgical gut suture is placed in tissue, a moderate tissue inflammation occurs characteristic of foreign body response to a substance. This is followed by a loss of tensile strength and suture mass, as the proteolytic enzymatic digestive process dissolves the surgical gut. This process continues until the suture is completely absorbed. Many variable factors may affect the rate of absorption. Some of the major factors which can affect tensile strength loss and absorption rates are:

1. Type of suture-Plain gut is expected to absorb more rapidly than chromic gut.
2. Infection-Surgical gut is absorbed more rapidly in infected tissue than in non-infected tissue.
3. Tissue Sites-Surgical gut will absorb more rapidly in tissue where increased levels of proteolytic enzymes are present, as in the secretions exhibited in the stomach, cervix and vagina.

CONTRAINDICATIONS

The use of this suture is contraindicated in patients with known sensitivities or allergies to collagen or chromium, as gut is a Collagen based material, and chromic gut is treated with chromic salt solutions.

WARNINGS

Prolonged contact of this or any other suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation. This suture, being absorbable, should not be used where extended approximation of tissue is required.

Under some circumstances, notably orthopedic procedures, immobilization by external support may be employed at the discretion of the surgeon. Do not resterilize. Discard open, unused sutures. Store at room temperature.

Avoid prolonged exposure to elevated temperatures.

Users should be familiar with surgical procedures and techniques involving gut suture before using Plain and Chromic Gut for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used.

Since any foreign material in the presence of bacterial contamination may enhance bacterial infectivity, acceptable surgical practice should be followed with respect to drainage and closure of contaminated or infected wounds.

The use of this suture may be inappropriate in elderly, malnourished, or debilitated patients, or in patients suffering from cancer, anemia, obesity, diabetes, infection, or other conditions which may delay wound healing.

As this is an absorbable material, the use of supplemental nonabsorbable sutures should be considered by the surgeon in the closure of the abdomen, chest, joints or other sites which may undergo expansion, stretching, or distention of which may require additional support.

Certain patients may be hypersensitive to collagen or chromium and might exhibit an immunological reaction resulting in inflammation, tissue granulation or fibrosis, wound suppuration and bleeding, as well as sinus formation.

PRECAUTIONS

In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

The surgeon should avoid unnecessary tension when running down knots, to reduce the occurrence of surface fraying and weakening of the strand.

ADVERSE REACTIONS

Adverse effects associated with the use of this device include: wound dehiscence, increased rates of absorption over time (depending on the type of suture used, the presence of infection and the tissue site), failure to provide adequate wound support in closure of the abdomen, chest joints, and other sites where expansion, stretching or distension occur, etc., unless addition support-is supplied through the use of nonabsorbable suture material, failure to provide adequate wound support in elderly, malnourished or debilitated patients or in patients suffering from cancer, anemia, obesity, diabetes, infection or other conditions which may delay wound healing, allergic response in patients with known sensitivities to collagen or chromium which may result in an immunological reaction resulting in inflammation, tissue granulation or fibrosis, wound suppuration and bleeding, as well as sinus formation, enhanced bacterial infectivity, moderate tissue inflammatory response characteristic of foreign body response, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, and pain, edema and erythema at the wound site.

HOW SUPPLIED

Plain and Chromic Gut is available in sizes 6-0 through 3-0. The sutures are available plain or chromic and packaged in a 89% isopropanol, 10% water, and 1% triethanolamine tubing fluid. The Chromic Gut has been tanned by treatment with chromium salts in trivalent form and with oxidised pyrogallol (oxidised pyrogallic acid) to prolong its resistance to absorption.

The suture is supplied sterile in pre-cut lengths or ligating reels, both non-needled or affixed to various needle types; one-dozen boxes.

Caution: National or local regulations may restrict this device to sale, distribution and use by or on the order of a physician, veterinarian or other licensed practitioner.

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