POLYPROPYLENE

Instructions for Use

POLYPROPYLENE (NONABSORBABLE SURGICAL SUTURES, USP)

DESCRIPTION

Polypropylene is a nonabsorbable, sterile surgical suture composed of a strand of polypropylene, a synthetic linear polyolefin. The suture is provided dyed blue.

The pigment for the blue dyed suture is phthalocyaninato (2-) Copper. Polypropylene meets all requirements established by the United States Pharmacopeia for Nonabsorbable Surgical Sutures.

INDICATIONS

Polypropylene is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neurological procedures.

ACTIONS

Polypropylene elicits a minimal acute inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissue. Polypropylene is not absorbed, nor is any significant change in strength retention known to occur in vivo.

CONTRAINDICATIONS

None known.

WARNINGS

Prolonged contact of this or any other suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation. Do not resterilize. Discard open, unused sutures.

Users should be familiar with surgical procedures and techniques involving nonabsorbable sutures before employing Polypropylene for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used.

As any foreign material in the presence of bacterial contamination may enhance bacterial infectivity, acceptable surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds.

PRECAUTIONS

In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

As with any suture material, adequate knot security requires the accepted surgical technique of flat, square ties with additional throws as warranted by surgical circumstance and the experience of the surgeon. The use of additional throws may be particularly appropriate when knotting monofilaments.

ADVERSE REACTIONS

Adverse effects, associated with the use of this device include: wound dehiscence, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, enhanced infectivity, minimal acute inflammatory tissue reaction, and pain, edema, and erythema at the wound site.

HOW SUPPLIED

Polypropylene is available in sizes 10-0 through 2. The suture is provided dyed blue.

The pigment for the blue dyed suture is phthalocyaninato (2-) Copper. The suture is supplied sterile in pre-cut lengths or ligating reels; both non-needled and affixed to various needle types; in one-dozen boxes.

Caution: National or local regulations may restrict this device to sale, distribution and use by or on the order of a physician, veterinarian or other licensed practitioner.

ProAdvantage by NDC
Made in BUSAN, KOREA

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