

Reusable

Hot Packs & Terry Covers

Reusable Hot Packs

- Designed to provide safe and effective moist heat for pain relief or therapy for sore, stiff muscles and joints
- Provides up to 30 minutes of soothing moist heat that can improve blood circulation, reduce inflammation and increase flexibility
- Reusable design for numerous effective and economical treatments for professional or home use
- Immerse in hot water or heating unit
- Use with All-Terry or Foam Filled terry covers
- Several sizes available to adapt to the contours of the body
- Constructed of a synthetic fabric and filled with bentonite

Covers

- Custom designed for reusable hot packs
- Choice of all-terry or foam filled terry covers
- · Washable and reusable




Hot Pack with Cover

Distributed by:



Reusable

Hot Packs & Terry Covers





Item #	Description	Item#	Description
--------	-------------	-------	-------------

REUSABLE HOT PACKS

 P503012
 Standard, 10" x 12"
 P503018
 Spinal, 10" x 18"

 P503015
 Oversized, 15" x 24"
 P503010
 Spinal, 10" x 24"

P503024 Cervical Contour, 24" long P503020 Knee & Shoulder with Hole, 10" x 20"







Item # Description Item # Description

TERRY COVERS

Foam Filled		All-Terry	All-Terry	
P503130	Standard, 10" x 28"	P503110	Standard, 10" x 28"	
P503140	Oversized, 24" x 32"	P503120	Oversized, 24" x 32"	
P503150	Cervical Contour, 9" x 26"	P503100	Cervical Contour, 9" x 26"	