

Urine

Reagent Strips

Pro
ADVANTAGE[®]

by **NDC**

Pro Advantage Urine Reagent Strips allow healthcare professionals to detect signs of infection, urinary tract damage, kidney and liver disease, fat breakdown, diabetes, and more in about a minute. Results are obtained by direct comparison of the test strip with the color blocks printed on the bottle label or with use on the Pro Advantage Urine Analyzer. Available in 10, 7, 4 and 2 parameters.

- **10 parameter Urine Reagent Strips**

For the semi-quantitative and qualitative detection of Leukocytes, Nitrite, Urobilinogen, Protein, pH, Blood, Specific Gravity, Ketone, Bilirubin, Glucose in urine.

- **7 parameter Urine Reagent Strips**

For the semi-quantitative and qualitative detection of Leukocytes, Nitrite, Protein, pH, Blood, Ketone, Glucose in urine.

- **4 parameter Urine Reagent Strips**

For the semi-quantitative and qualitative detection of Leukocytes, Nitrite, Protein, Glucose in urine.

- **2 parameter Urine Reagent Strips**

For the semi-quantitative and qualitative detection of glucose and protein in urine.



Item #	Description	Bottle
P080010	10 parameter Urine Reagent Strips	100
P080017	7 parameter Urine Reagent Strips	100
P080014	4 parameter Urine Reagent Strips	100
P080012	2 parameter Urine Reagent Strips	100



For Urinalysis, Visual
CPT Code: 81002
For Urinalysis, Automated
CPT Code: 81003QW

Distributed By:

Urine

Analyzer

Pro
ADVANTAGE[®]

by **NDC**

Portable urine chemistry analyzer exclusively for reading Pro Advantage Urine Reagent Strips.

- Process 120 tests/hour
- Large 2,000 patient memory storage
- EMR/LIS Connectivity
- Eliminates color interpretation variability
- Minimizes transcription errors with an easy-to-read printout
- Patient ID Capable
- LCD screen
- User friendly, portable and convenient



Item #	Description
P080000	Urine Analyzer
P080001	Barcode Reader
P080002	Liquid Controls <i>Contains: 1-5 mL Bottle Positive Control, 1-5 mL Bottle Negative Control</i>
P080004	Thermal Printer Paper, 4 roll/bx
P080005	Data Transfer Cable



For Urinalysis, Visual
CPT Code: 81002
For Urinalysis, Automated
CPT Code: 81003QW

Distributed By: